

# Silicone Gastrostomy Tube



## Indications

This product is suitable for the replacement of the existing device of the established stoma tube, and is placed in the stomach to realize the delivery of enteral nutrient solution; gastric decompression.



## Specifications

Tube Specification	Rated Volume of Balloon (ml)	Tube Length (cm)
12 Fr	5	0.8/1.0/1.2/1.5 1.7/2.0/2.3/2.5 2.7/3.0/3.5/4.0 4.5/5.0/5.5/6.0 6.5
14 Fr	5	
16 Fr	10	
18 Fr	10	
20 Fr	10	
22 Fr	10	
24 Fr	10	



## Specifications

Tube Specification	Rated Volume of Balloon (ml)	Tube Length (mm)	Mark Location (cm)	Color Code
12 Fr	5	155	2 / 4 / 6	white
14 Fr	5			green
16 Fr	20			orange
18 Fr	20			red
20 Fr	20			yellow
22 Fr	20			purple
24 Fr	20			navyblue