

Medical Personal Protection System

Indications

It can be worn by medical staff in operating rooms, wards and diagnosis and treatment areas to help prevent the spread of infection sources to patients and the spread of infection sources carried in patients' blood or body fluids to medical staff.

Features

- ※ The dual airflow system improves the cooling effect and can realize the adjustment of the airflow from 170L/min to 500L/min.
- ※ The battery provides up to 10+ hours of working time.
- ※ All-in-one full body protection, safer breathable barrier to body fluids and viruses, meeting AAMI level 4 and EN13795 requirements.
- ※ Wide-angle anti-fog screen, clear visor, easy to operate.
- ※ Lighter materials and a well-balanced center of gravity reduce muscle fatigue.



Accessories

- ※ Rechargeable Battery
- ※ 4-Station Charger
- ※ Lightweight LED light
- ※ Disposable sterile headgear

Components

This product consists of helmet, wall rack, power pack, battery charger, hood and gown.

Specifications

HP 18A